

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Spin Echo
Coil Configuration	HD 8Ch Body Array by GE	Imaging Options	Seq, Fast, SS
Plane	3-PLANE	SCAN RANGE	
Series Description	3 PLANE LOC	FOV	42.0
SCAN TIMING		Slice Thickness	8.0
TE	80.0	Slice Spacing	0.0
TR	Minimum	ACQ TIMING	
Receiver Bandwidth	83.33	Freq	320
IMAGE ENHANCE		Phase	160
Filter Choice	None	Freq DIR	Unswap
GATING/TRIGGER		# of Acq. Before Pause	0
Auto Trigger Type	Off	Phase FOV	1.00
FMRI		Auto Shim	Auto
PSD Trigger	Internal	Phase Correction	No
Slice Order	Interleaved	USER CVS	
View Order	Bottom/Up	User CV1	1.00
# of Repetitions REST	0	User CV2	240.00
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		Seperate Series	0
Tag Type	None	Mask Phase	0
TRICKS		Mask Pause	0
Pause On/Off	On	DIFFUSION	
Auto Subtract	0	Recon All Images	On
Auto SCIC	Off	CONTRAST	
		Contrast Yes/No	No

3 PLANE LOC

3 PLANE LOC

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	<i>Feet First</i>	Imaging Mode	<i>3D</i>
Patient Position	<i>Supine</i>	Pulse Sequence	<i>SPGR</i>
Coil Configuration	<i>HD 8Ch Body Array by GE</i>	Imaging Options	<i>EDR, Fast, ZIP2, Calib</i>
Plane	<i>AXIAL</i>	IMAGE ENHANCE	
Series Description	<i>Calibration Scan</i>	Filter Choice	<i>None</i>
SCAN RANGE		GATING/TRIGGER	
FOV	<i>50.0</i>	Auto Trigger Type	<i>Off</i>
Slice Thickness	<i>16.0</i>	MULTI-PHASE	
ACQ TIMING		Seperate Series	<i>0</i>
Freq DIR	<i>R/L</i>	Mask Phase	<i>0</i>
Auto Shim	<i>Auto</i>	Mask Pause	<i>0</i>
Phase Correction	<i>No</i>	DIFFUSION	
FMRI		Recon All Images	<i>On</i>
PSD Trigger	<i>Internal</i>	CONTRAST	
Slice Order	<i>Interleaved</i>	Contrast Yes/No	<i>No</i>
View Order	<i>Bottom/Up</i>		
# of Repetitions REST	<i>0</i>		
# of Repetitions ACTIVE	<i>0</i>		
SAT			
Tag Type	<i>None</i>		
TRICKS			
Pause On/Off	<i>On</i>		
Auto Subtract	<i>0</i>		
Auto SCIC	<i>Off</i>		

Calibration Scan

Calibration Scan

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Spin Echo
Coil Configuration	HD 8Ch Body Array by GE	Imaging Options	EDR, Fast, SS, Asset
Plane	AXIAL	SCAN RANGE	
Series Description	AX T2 (SSFSE)	FOV	38.0
SCAN TIMING		Slice Thickness	5.0
TE	140.0	Slice Spacing	1.0
TR	1000.0	ACQ TIMING	
Receiver Bandwidth	83.33	Freq	320
IMAGE ENHANCE		Phase	224
Filter Choice	B	Freq DIR	R/L
GATING/TRIGGER		# of Acq. Before Pause	22
Auto Trigger Type	Off	Phase FOV	0.90
FMRI		Auto Shim	Auto
PSD Trigger	Internal	Phase Correction	No
Slice Order	Interleaved	USER CVS	
View Order	Bottom/Up	User CV2	240.00
# of Repetitions REST	0	User CV15	1.00
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		Seperate Series	0
Tag Type	None	Mask Phase	0
TRICKS		Mask Pause	0
Pause On/Off	On	DIFFUSION	
Auto Subtract	0	Recon All Images	On
Auto SCIC	2	CONTRAST	
		Contrast Yes/No	No

AX T2 (SSFSE)

AX T2 (SSFSE)

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Spin Echo
Coil Configuration	HD 8Ch Body Array by GE	Imaging Options	EDR, Fast, SS, Asset
Plane	CORONAL	SCAN RANGE	
Series Description	COR T2 (SSFSE)	FOV	38.0
SCAN TIMING		Slice Thickness	5.0
TE	140.0	Slice Spacing	1.0
TR	1000.0	ACQ TIMING	
Receiver Bandwidth	83.33	Freq	320
IMAGE ENHANCE		Phase	224
Filter Choice	B	Freq DIR	S/I
GATING/TRIGGER		# of Acq. Before Pause	18
Auto Trigger Type	Off	Phase FOV	1.00
FMRI		Auto Shim	Auto
PSD Trigger	Internal	Phase Correction	No
Slice Order	Interleaved	USER CVS	
View Order	Bottom/Up	User CV2	240.00
# of Repetitions REST	0	User CV15	1.00
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		Seperate Series	0
Tag Type	None	Mask Phase	0
TRICKS		Mask Pause	0
Pause On/Off	On	DIFFUSION	
Auto Subtract	0	Recon All Images	On
Auto SCIC	2	CONTRAST	
		Contrast Yes/No	No

COR T2 (SSFSE)

COR T2 (SSFSE)

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	FRFSE-XL
Coil Configuration	HD 8Ch Body Array by GE	Imaging Options	FC, NPW, EDR, TRF, Fast, ZIP512, FR
Plane	AXIAL	IDEAL	3
Series Description	AX T2 (SMALL FOV)		
SCAN TIMING		SCAN RANGE	
Flip Angle	160	FOV	20.0
TE	110.0	Slice Thickness	5.0
Number of Echoes	1	Slice Spacing	1.0
TR	3989.0		
Echo Train Length	18	ACQ TIMING	
Receiver Bandwidth	31.25	Freq	320
		Phase	256
IMAGE ENHANCE		Freq DIR	R/L
Filter Choice	B	Fat Shift DIR	Normal (L)
GATING/TRIGGER		NEX	3.00
Auto Trigger Type	Off	Auto Shim	Auto
FMRI		Phase Correction	No
PSD Trigger	Internal	Flow Direction Compensation	Slice
Slice Order	Interleaved	USER CVS	
View Order	Bottom/Up	TR Min	2500.0
# of Repetitions REST	0	TR Max	11000.0
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		Seperate Series	0
SAT Location	A	Mask Phase	0
SAT Location	S	Mask Pause	0
SAT Location	I	DIFFUSION	
Tag Type	None	Recon All Images	On
TRICKS		CONTRAST	
Pause On/Off	On	Contrast Yes/No	No
Auto Subtract	0		
Auto SCIC	2		


AX T2 (SMALL FOV)

AX T2 (SMALL FOV)

OTHERS

Protocol Notes

STRAIGHT AXIAL
S/I/A SAT BANDS LIKE
BELOW
FOR UTERINE ANOMALY
(Angle slices parallel to
uterine strip like below)



BH AX 3D LAVA (PELVIS)	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>3D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>LAVA-Flex</i>
	Coil Configuration	<i>HD 8Ch Body Array by GE</i>	Imaging Options	<i>EDR, Fast, ZIP512, ZIP2, ARC</i>
	Plane	<i>OBLIQUE</i>	IDEAL	<i>7</i>
	Series Description	<i>BH AX 3D LAVA (PELVIS)</i>	SCAN RANGE	
	SCAN TIMING		FOV	<i>38.0</i>
	Flip Angle	<i>12</i>	Slice Thickness	<i>4.0</i>
	TE	<i>Min Full</i>	Location per Slab	<i>48</i>
	Receiver Bandwidth	<i>83.33</i>	Overlap Locations	<i>0</i>
	IMAGE ENHANCE		ACQ TIMING	
	Filter Choice	<i>B</i>	Freq	<i>288</i>
	GATING/TRIGGER		Phase	<i>224</i>
	Auto Trigger Type	<i>Off</i>	Freq DIR	<i>R/L</i>
	MULTI-PHASE		NEX	<i>1.00</i>
	Seperate Series	<i>0</i>	Phase FOV	<i>0.80</i>
	Delay after Acquisition without AV	<i>20</i>	Auto Shim	<i>Auto</i>
	Trigger Delay without AV	<i>0</i>	Phase Correction	<i>No</i>
	Mask Phase	<i>0</i>	FMRI	
	Mask Pause	<i>0</i>	PSD Trigger	<i>Internal</i>
	DIFFUSION		Slice Order	<i>Interleaved</i>
	Recon All Images	<i>On</i>	View Order	<i>Bottom/Up</i>
	CONTRAST		# of Repetitions REST	<i>0</i>
	Contrast Yes/No	<i>No</i>	# of Repetitions ACTIVE	<i>0</i>
			SAT	
		Tag Type	<i>None</i>	
		TRICKS		
		Pause On/Off	<i>On</i>	
		Auto Subtract	<i>0</i>	
		Auto SCIC	<i>2</i>	
		OTHERS		
		Protocol Notes	<i>Breath hold Straight Axial covering entire pelvis</i>	

BH AX 3D LAVA (PELVIS)

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Spin Echo
Coil Configuration	HD 8Ch Body Array by GE	Imaging Options	EDR, Fast, SS, Asset
Plane	SAGITTAL	SCAN RANGE	
Series Description	SAG T2 (REST)	FOV	30.0
SCAN TIMING		Slice Thickness	6.0
TE	140.0	Slice Spacing	1.0
TR	1000.0	ACQ TIMING	
Receiver Bandwidth	83.33	Freq	320
IMAGE ENHANCE		Phase	224
Filter Choice	B	Freq DIR	S/I
GATING/TRIGGER		# of Acq. Before Pause	0
Auto Trigger Type	Off	Phase FOV	0.80
FMRI		Auto Shim	Auto
PSD Trigger	Internal	Phase Correction	No
Slice Order	Interleaved	USER CVS	
View Order	Bottom/Up	User CV2	240.00
# of Repetitions REST	0	User CV15	1.00
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		Seperate Series	0
Tag Type	None	Mask Phase	0
TRICKS		Mask Pause	0
Pause On/Off	On	DIFFUSION	
Auto Subtract	0	Recon All Images	On
Auto SCIC	2	CONTRAST	
OTHERS		Contrast Yes/No	No
Protocol Notes	LRG FOV TO INCLUDE KIDNEYS		

SAG T2 (REST)

SAG T2 (REST)

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Spin Echo
Coil Configuration	HD 8Ch Body Array by GE	Imaging Options	EDR, Fast, SS, Asset
Plane	SAGITTAL	SCAN RANGE	
Series Description	SAG T2 (STRAIN)	FOV	30.0
SCAN TIMING		Slice Thickness	6.0
TE	140.0	Slice Spacing	1.0
TR	1000.0	ACQ TIMING	
Receiver Bandwidth	83.33	Freq	320
IMAGE ENHANCE		Phase	224
Filter Choice	B	Freq DIR	S/I
GATING/TRIGGER		# of Acq. Before Pause	0
Auto Trigger Type	Off	Phase FOV	0.80
FMRI		Auto Shim	Auto
PSD Trigger	Internal	Phase Correction	No
Slice Order	Interleaved	USER CVS	
View Order	Bottom/Up	User CV2	240.00
# of Repetitions REST	0	User CV15	1.00
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		Seperate Series	0
Tag Type	None	Mask Phase	0
TRICKS		Mask Pause	0
Pause On/Off	On	DIFFUSION	
Auto Subtract	0	Recon All Images	On
Auto SCIC	2	CONTRAST	
OTHERS		Contrast Yes/No	No
Protocol Notes	ASK PT TO BEAR DOWN TO STRAIN WHEN SCAN BEGINS		

SAG T2 (STRAIN)

SAG T2 (STRAIN)

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Spin Echo
Coil Configuration	HD 8Ch Body Array by GE	Imaging Options	EDR, Fast, SS, Asset
Plane	SAGITTAL	SCAN RANGE	
Series Description	SAG T2 (SQUEEZE)	FOV	30.0
SCAN TIMING		Slice Thickness	6.0
TE	140.0	Slice Spacing	1.0
TR	1000.0	ACQ TIMING	
Receiver Bandwidth	83.33	Freq	320
IMAGE ENHANCE		Phase	224
Filter Choice	B	Freq DIR	S/I
GATING/TRIGGER		# of Acq. Before Pause	0
Auto Trigger Type	Off	Phase FOV	0.80
FMRI		Auto Shim	Auto
PSD Trigger	Internal	Phase Correction	No
Slice Order	Interleaved	USER CVS	
View Order	Bottom/Up	User CV2	240.00
# of Repetitions REST	0	User CV15	1.00
# of Repetitions ACTIVE	0	MULTI-PHASE	
SAT		Seperate Series	0
Tag Type	None	Mask Phase	0
TRICKS		Mask Pause	0
Pause On/Off	On	DIFFUSION	
Auto Subtract	0	Recon All Images	On
Auto SCIC	2	CONTRAST	
OTHERS		Contrast Yes/No	No
Protocol Notes	ASK PT TO SQUEEZE RECTUM WHEN SCAN BEGINS		

SAG T2 (SQUEEZE)

SAG T2 (SQUEEZE)

PATIENT POSITION		IMAGING PARAMETERS	
Patient Entry	Feet First	Imaging Mode	2D
Patient Position	Supine	Pulse Sequence	Fiesta
Coil Configuration	HD 8Ch Body Array by GE	Imaging Options	Seq, EDR, Fast, MPH
Plane	SAGITTAL	SCAN RANGE	
Series Description	SAG FIESTA (STRAIN)	FOV	28.0
SCAN TIMING		Slice Thickness	8.0
Flip Angle	55	Slice Spacing	0.0
TE	Min Full	ACQ TIMING	
Receiver Bandwidth	83.33	Freq	140
IMAGE ENHANCE		Phase	192
Filter Choice	None	Freq DIR	S/I
GATING/TRIGGER		NEX	1.00
Auto Trigger Type	Off	Phase FOV	0.90
FMRI		Auto Shim	Auto
PSD Trigger	Internal	Phase Correction	No
Slice Order	Interleaved	USER CVS	
View Order	Bottom/Up	User CV20	2.00
# of Repetitions REST	0	MULTI-PHASE	
# of Repetitions ACTIVE	0	Slice per Location	30
DIFFUSION		Seperate Series	0
Recon All Images	On	Mask Phase	0
CONTRAST		Mask Pause	0
Contrast Yes/No	No	TRICKS	
		Pause On/Off	On
		Auto Subtract	0
		Auto SCIC	Off
		OTHERS	
		Protocol Notes	ONE SLICE MIDLINE ASK PT TO BEAR DOWN TO STRAIN WHEN SCAN BEGINS

SAG FIESTA (STRAIN)

SAG FIESTA (STRAIN)

Protocol: adult_pelvis_DEFECOGRAPHY_20210521081344861_1

SAG FIESTA (EVACUATE)	PATIENT POSITION		IMAGING PARAMETERS	
	Patient Entry	<i>Feet First</i>	Imaging Mode	<i>2D</i>
	Patient Position	<i>Supine</i>	Pulse Sequence	<i>Fiesta</i>
	Coil Configuration	<i>HD 8Ch Body Array by GE</i>	Imaging Options	<i>Seq, EDR, Fast, MPh</i>
	Plane	<i>SAGITTAL</i>	SCAN RANGE	
	Series Description	<i>SAG FIESTA (EVACUATE)</i>	FOV	<i>28.0</i>
	SCAN TIMING		Slice Thickness	<i>8.0</i>
	Flip Angle	<i>55</i>	Slice Spacing	<i>0.0</i>
	TE	<i>Min Full</i>	ACQ TIMING	
	Receiver Bandwidth	<i>83.33</i>	Freq	<i>140</i>
	IMAGE ENHANCE		Phase	<i>192</i>
	Filter Choice	<i>None</i>	Freq DIR	<i>S/I</i>
	GATING/TRIGGER		NEX	<i>1.00</i>
	Auto Trigger Type	<i>Off</i>	Phase FOV	<i>0.90</i>
	FMRI		Auto Shim	<i>Auto</i>
	PSD Trigger	<i>Internal</i>	Phase Correction	<i>No</i>
	Slice Order	<i>Interleaved</i>	USER CVS	
	View Order	<i>Bottom/Up</i>	User CV20	<i>2.00</i>
	# of Repetitions REST	<i>0</i>	MULTI-PHASE	
	# of Repetitions ACTIVE	<i>0</i>	Slice per Location	<i>30</i>
	DIFFUSION		Separate Series	<i>0</i>
	Recon All Images	<i>On</i>	Mask Phase	<i>0</i>
	CONTRAST		Mask Pause	<i>0</i>
	Contrast Yes/No	<i>No</i>	TRICKS	
			Pause On/Off	<i>On</i>
		Auto Subtract	<i>0</i>	
		Auto SCIC	<i>Off</i>	
		OTHERS		
		Protocol Notes	<i>ONE SLICE MIDLINE</i>	

SAG FIESTA (EVACUATE)